

SS-Ethylendiaminodisuccinic Acid (EDDS)

0.5 mg/ml (H₂O)

HP1090 Liquid Chromatograph with diode array detector

Column: Equisil C-18 (5 µm) 125 + 20 mm x 3 mm I.D.

Eluent A: 20 mM Na/K-phosphate buffer pH 7.2 + 5 mM tetrabutylammonium-hydrogensulfate (TBA)

Eluent B: acetonitrile

Gradient profile:

0 min	4% B
5 min	4% B
15 min	60% B
16 min	60% B

Flow rate: 0.85 ml/min

Injection volume: 10 µl

